

Monoclonal antibodies Aerosol formulations

Renaud Respaud - 02/07/2015



MAbDelivery
Industrial Workshop
www.mabdelivery.fr

July 2nd 2015, Tours
VINCI International Convention Centre - France

POLEPHARMA **MAbiimprove**
LabEx Tours Montpellier


Challenges for mAbs inhalation therapy

Effect of formulation on the stability of a nebulized antibody

Technical challenges for a successful inhalation therapy

Dependent on:

- Aerosol technology
- Performance of the device (aerosol output, particle size)



Particle size



> 5 μ m Upper
Respiratory tract

2- 5 μ m Mid Respiratory tract

Respirable
0.5-2 μ m Alveoli

Technical challenges for a successful inhalation therapy



Dependent on:

- Aerosol technology
- Performance of the device (aerosol output, particle size)
- Drugs
- Characteristics of drug formulation (concentration, pH, viscosity...)


Technical challenges for a successful inhalation therapy



Dependent on:

- Aerosol technology
- Performance of the device (aerosol output, particle size)
- Drugs
- Characteristics of drug formulation (concentration, pH, viscosity...)
- Loss of drug along the aerosolisation process

Technical challenges for a successful inhalation therapy



Dependent on:

- Aerosol technology
- Performance of the device (aerosol output, particle size)
- Drugs
- Characteristics of drug formulation (concentration, pH, viscosity...)
- Loss of drug along the aerosolisation process

Strategy of development for a formulation for mAbs inhalation

1. Solid or liquid aerosols ?

Solid or liquid aerosols ?

Solid aerosols: DPIs - Dry Powder Inhalers

- **Advantages** : user friendly / no cleaning / quick delivery
- **Disadvantages** : manufacturing process (gas-liquid interface, thermal stress) / drug and device development / Quantity of drugs



Liquid aerosols:

- MDIs - Metered Dose Inhalers
- Nebulizers
 - US
 - Jet
 - Mesh



Strategy of development for a formulation for mAbs inhalation

1. Solid or liquid aerosols 

2. Type of nebulizer ?

Type of nebulizer ?







Figure 2. Mechanisms of action of nebulizers.

Strategy of development for a formulation for mAbs inhalation

1. Solid or liquid aerosols 
2. Type of nebulizer 
3. Limitation of nebulizer regarding mAbs ?

Limitation of nebulizer = physical stresses

Mesh nebulizer



Dispersion/suspension of liquid droplets in a gaseous medium

- Heat
- Shearing
- Air-liquid interface ++++++
- +

Respaud *et al.* 2015

Strategy of development for a formulation for mAbs inhalation


1. Solid or liquid aerosols 

2. Type of nebulizer 

3. Limitation of nebulizer regarding mAbs 

4. Formulation challenges ?

Formulation Challenges







Prone to various types of physical and chemical degradation

- Denaturation, deamidation, oxidation
- **Non-covalent/covalent aggregation +++++**

- **Antibody concentration:** known to affect aggregation
- **Buffer dilution and pH:** PBS, citric acid, histidine...
- **Excipients:**
 - **Surfactants (Polysorbate...):** protect mAbs at the air-liquid interface
 - **Sugars (sucrose, trehalose...):** cryoprotection
 - **Amino-acids (glycine, lysine, isoleucine...):** stabilizers

Strategy of development for a formulation for mAbs inhalation

1. Solid or liquid aerosols 
2. Type of nebulizer 
3. Limitation of nebulizer regarding mAbs 
4. Formulation challenges 


Challenges for mAbs inhalation therapy

Effect of formulation on the stability of a nebulized antibody

Context

- Drug & device for the rapid delivery into the alveoli of anti-ricin mAb
 - Drug : IgG1 43RCA-G1 in PBS pH 7.2
 - Device : mesh nebulizer
 - Formulation ???
- Aggregation of a nebulized polyclonal IgG and IgG1 43RCA-G1
 - Effect of Antibody concentration
 - Influence of the device
 - Effect of surfactants

High concentration seems to protect polyclonal antibodies



Fluorescence Microscopy (Nile Red)
Large aggregates > 2 μ m



DLS: Dynamic Light Scattering
Medium-size aggregates < 1 μ m

Antibody concentration (mg/ml)	Hydrodynamic diameter (nm \pm SD)	Percentage monomers (mean% mass)
2	NA	NA
5	NA	NA
10	NA	NA
15	12.6 \pm 1.42	94.4 \pm 9.0
30	11.9 \pm 0.06	99.6 \pm 0.2
40/45	12.6 \pm 0.35	99.6 \pm 0.5

Device has influence on large aggregates

Table 1. Characteristics of a 10 mg/ml solution of antibody in PBS pH 7.2 after nebulization with three nebulizers. Each result is the mean (\pm SD) of three nebulizations

Nebulizer	Hydrodynamic diameter (nm \pm SD)	Percentage monomers (% mass)	Size exclusion chromatography (% monomers/% other aggregates)	VMD (μ m) (mean \pm SD)	Flow rate (ml/min \pm SD)
Omron	NA	NA	91.1 \pm 0.6 / 8.9 \pm 0.5	5.6 \pm 0.1	0.23 \pm 0.04
Aerogen Solo	NA	NA	91.3 \pm 2.5 / 8.7 \pm 1.9	4.1 \pm 0.1	0.50 \pm 0.01
PARI eFlow	NA	NA	90.6 \pm 0.5 / 9.5 \pm 0.6	3.8 \pm 0.3	0.45 \pm 0.05




Respaud R et al. 2014

Surfactants protect polyclonal antibodies during nebulization

PS-20 concentration% (w/v)	Hydrodynamic diameter (nm ± SD)	Percentage monomers (% mass ± SD)
0%	NA	NA
0.00001%	NA	NA
0.0001%	NA	NA
0.001%	11.4 ± 0.5	99.9 ± 0.2
0.01%	11.2 ± 0.3	99.9 ± 0.2

DLS: Dynamic Light Scattering
Medium-size aggregates < 1µm




Fluorescence Microscopy (Nile Red)
Large aggregates > 2µm

Does surfactants protect monoclonal antibodies during nebulization ?

Characteristics of solutions of IgG 43RCA-G1 15 mg/mL in buffer 25 mM, NaCl 135 mM, with various concentration of PS20, before and after nebulization. Each result is the result of 1 nebulization


Nebulizer	Surfactant (PS20) concentration (% w/v)	SEC (% of monomer / % of other species)		DLS (% of monomer)		Large antibody aggregates (10 ³ /mL)	
		Before	After	Before	After	Before	After
		nebulization	nebulization	nebulization	nebulization	nebulization	nebulization
Aerogen	0	100 / 0	100 / 0	99.6	99.1	85.6	128.8
solo	0.01	100 / 0	99.9 / 0.1	99.2	98.6	27.9	45.9
	0.1	100 / 0	99.8 / 0.2	99.9	100	33.3	29.7

No dimer (<0.5%) of IgG1 43RCA-G1 during nebulization



Size exclusion chromatography of a 15 mg/ml 43RCA-G1 IgG1 solution in developed buffer with various concentration of PS20 before and after nebulization with the Aerogen Solo. UV signal

Does surfactants protect monoclonal antibodies during nebulization ?



Drug & device for the rapid delivery into the alveoli of anti-ricin mAb

Conclusion and perspectives

- Utility of **adding surfactants** and **increasing protein concentration** to stabilize antibody formulations during nebulization
- Optimization for each **“drug and device” pairing** (specifications of the device and **properties of the mAb molecule**)
- Chemical modifications of mAbs during nebulization
- Toxicity of excipients used to stabilize mAbs during nebulization

Thank you for your attention



Plateforme de physicochimie
des Biomédicaments





Nathalie Heuzé-Vourc'h
Laurent Vecellio
Denis Marchand
Christelle Parent
Jeffrey Pardessus
Patrice Diot
Agnès Maillet



Hervé Watier
André Pèlerin
Marc Bonnemaïson
Annie Gauvineau
Arnaud Daguet
Alexandra Farrell





Model of the spatial distribution of antibodies in aerosol droplets as a function of antibody concentration

